

# The Education University of Hong Kong

## 2021-2022 Quality Education Fund Thematic Network – Tertiary Institutes



### STEM Project Team

SCHOOL: CHRISTIAN AND MISSIONARY ALLIANCE SUN  
KEI SECONDARY SCHOOL (S4)



TOPIC: WORKSHOP 1 – DEVELOP PORTABLE  
SPECTROPHOTOMETER WITH DATA LOGGING

# REQUIREMENT LIST:

- - SEN-1 5050 SPARKFUN TRIAD SPECTROSCOPY SENSOR -AS7265X
- - ARDUINO-UNO (WITH POWER SUPPLY)
- - COMPUTER WITH ARDUINO IDE AND EXCEL
- - SAMPLE HOLDER
- - MIRROR
- - RUBBER BANDS
- - BOX(FOR COVERING THE WHOLE DEVICE FROM ANY OTHER LIGHT SOURCE)
- - SAMPLE WATER FOR DEMONSTRATION(RECOMMEND RELATED TO TEACHING PLAN)



# CIRCUIT



# CODING



```
sketch_dec07a | Arduino 1.8.3
File Edit Sketch Tools Help
sketch_dec07a
void setup() {
  // put your setup code here, to run once:
}

void loop() {
  // put your main code here, to run repeatedly:
}

2 Arduino/Genuino Uno on COM3
```



# BASIC READING WITH LEDS

```
Example2_BasicReadingsWithLEDS
//
#include "SparkFun_AS7265X.h" //Click here to get the library: https://librarymanager/All#SparkFun\_AS7265X
AS7265X sensor;

#include <Wire.h>

void setup()
{
  Serial.begin(115200);
  Serial.println("AS7265x Spectral Triad Example");

  Serial.println("Point the Triad away and press a key to begin with illumination...");
  while (Serial.available() == false)
  {
    //Do nothing while we wait for user to press a key
  }
  Serial.read(); //Throw away the user's button

  if (sensor.begin() == false)
  {
    Serial.println("Sensor does not appear to be connected. Please check wiring. Freezing...");
    while (1)
    {
      ;
    }
  }

  sensor.disableIndicator(); //Turn off the hide status LED

  Serial.println("A, B, C, D, E, F, G, H, R, I, S, J, T, U, V, W, X, L");
}



void loop()
{
  sensor.takeMeasurementsWithSub(); //This is a hard wait while all 16 channels are measured

  Serial.print(sensor.getCalibratedA()); //410nm
  Serial.print(", ");
  Serial.print(sensor.getCalibratedB()); //435nm
  Serial.print(", ");
  Serial.print(sensor.getCalibratedC()); //460nm
  Serial.print(", ");
  Serial.print(sensor.getCalibratedD()); //485nm
  Serial.print(", ");
  Serial.print(sensor.getCalibratedE()); //510nm
  Serial.print(sensor.getCalibratedF()); //535nm
  Serial.print(", ");
  Serial.print(sensor.getCalibratedG()); //560nm
  Serial.print(", ");
  Serial.print(sensor.getCalibratedH()); //585nm
  Serial.print(", ");
  Serial.print(sensor.getCalibratedI()); //610nm
  Serial.print(", ");
  Serial.print(sensor.getCalibratedJ()); //645nm
  Serial.print(", ");
  Serial.print(sensor.getCalibratedK()); //680nm
  Serial.print(sensor.getCalibratedL()); //705nm
  Serial.print(", ");
  Serial.print(sensor.getCalibratedM()); //730nm
  Serial.print(", ");
  Serial.print(sensor.getCalibratedN()); //760nm
  Serial.print(", ");
  Serial.print(sensor.getCalibratedO()); //810nm
  Serial.print(sensor.getCalibratedP()); //860nm
  Serial.print(", ");
  Serial.print(sensor.getCalibratedQ()); //900nm
  Serial.print(sensor.getCalibratedR()); //940nm
  Serial.print(", ");
  Serial.println();
}
```


```
Example1_BasicReadingsWithLEDS
void loop()
{
  sensor.takeMeasurementsWithSub(); //This is a hard wait while all 16 channels are measured

  Serial.print(sensor.getCalibratedA()); //410nm
  Serial.print(", ");
  Serial.print(sensor.getCalibratedB()); //435nm
  Serial.print(", ");
  Serial.print(sensor.getCalibratedC()); //460nm
  Serial.print(", ");
  Serial.print(sensor.getCalibratedD()); //485nm
  Serial.print(", ");
  Serial.print(sensor.getCalibratedE()); //510nm
  Serial.print(", ");
  Serial.print(sensor.getCalibratedF()); //535nm
  Serial.print(", ");
  Serial.print(sensor.getCalibratedG()); //560nm
  Serial.print(", ");
  Serial.print(sensor.getCalibratedH()); //585nm
  Serial.print(", ");
  Serial.print(sensor.getCalibratedI()); //610nm
  Serial.print(", ");
  Serial.print(sensor.getCalibratedJ()); //645nm
  Serial.print(", ");
  Serial.print(sensor.getCalibratedK()); //680nm
  Serial.print(sensor.getCalibratedL()); //705nm
  Serial.print(", ");
  Serial.print(sensor.getCalibratedM()); //730nm
  Serial.print(", ");
  Serial.print(sensor.getCalibratedN()); //760nm
  Serial.print(", ");
  Serial.print(sensor.getCalibratedO()); //810nm
  Serial.print(sensor.getCalibratedP()); //860nm
  Serial.print(", ");
  Serial.print(sensor.getCalibratedQ()); //900nm
  Serial.print(sensor.getCalibratedR()); //940nm
  Serial.print(", ");
  Serial.println();
}
```

# TOOLS > MANAGE LIBRARIES (CTRL+SHIFT+M) LIBRARY – AS7265X



# DOWNLOAD AND TEST WITH SERIAL MONITOR (CTRL+SHIFT+M)




```
COM20
10:53:19.171 -> 0,2.98,20.73,19.41,4.75,1.45,1.02,0.48,3.29,0.84,2.06,0.41,0.00,0.00,0.85,1.12,0.00,AS7265x Spectral Tr: ^
10:53:19.312 -> A,B,C,D,E,F,G,H,R,I,S,J,T,U,V,W,K,L
10:53:19.640 -> 0.00,2.98,20.73,19.41,4.75,1.45,1.02,0.48,3.29,0.84,2.06,0.41,0.00,0.00,0.85,1.12,0.00,0.00,
10:53:20.484 -> 0.00,2.98,20.73,19.41,4.75,1.45,1.02,0.48,3.29,0.84,2.06,0.41,0.00,0.00,0.85,1.12,0.00,0.00,
10:53:21.328 -> 0.00,2.98,20.73,19.41,4.75,1.45,1.02,0.48,3.29,0.84,2.06,0.41,0.00,0.00,0.85,1.12,0.00,0.00,
10:53:22.125 -> 0.00,2.98,20.73,19.41,4.75,1.45,1.02,0.97,3.29,1.26,2.06,0.41,0.00,0.00,0.85,1.12,0.00,0.00,
10:53:22.968 -> 0.00,2.98,20.73,19.41,3.96,
```

Autoscroll  Show timestamp Newline 115200 baud Clear output

ager/A11#SparkFun\_AS7265X

# BEFORE DATA COLLECTION

**(1) Warm-up(5mins)  
before using\*\***



**(2) Clean before using**






# DURING DATA COLLECTION


## (3) Mark time and action

	A	B	C	D	
1	Time(From)	Time(To)	Event	Remark	
2	15:44	15:51	Warm up		
3	15:52	15:55	H2O	15:54:30	
4	15:55	15:57:36	Sample 7	15:58:30	
5	16:02:20	16:05	Sample 5	16:03:30	
6	16:08:36	16:12:40	Sample 3	16:10:30	
7	16:15:20	16:18:20	Sample 6	16:16:30	
8	16:20:41	16:23:41	Sample 4	16:22:30	
9	16:25:40	16:28:40	Sample 2	16:27:30	
10	16:29:00	16:32:00	ACR	16:31:30	
11	16:33:17	16:36:00	H2O (2)	16:34:30	
12					
13					

# SET-UP




# DATA COLLECTION



# AT LEAST TAKE 3 SET OF SAMPLE DATA (0%,L%,H%)

e.g. Water(0%)



# AFTER DATA COLLECTION TAKE DATA (COPY & PASTE)

(4) Copy all data  
(Copy & paste)

